

Candidate: _____ Date: _____

ID#: _____

Skill Drill 24-3 One-rescuer Bag-mask Ventilation on a Child

Task: Perform one-person bag-mask ventilation on a child.					
Performance Observations: The candidate shall be able to correctly perform one-person bag-mask ventilation on a child.					
Candidate Directive: "Properly perform one-person bag-mask ventilation on a child."					
No.	Task Steps	First Test		Retest	
		P	F	P	F
1.	Open the airway and insert the appropriate airway adjunct.				
2.	Hold the mask on the patient's face with a one-handed head tilt-chin lift technique (EC grip). Ensure a good mask-to-face seal while maintaining the airway.				
3.	Squeeze the bag 12 to 20 times a minute for infants and children. Ventilations should be delivered over 1 second. Allow adequate time for exhalation.				
4.	Assess effectiveness of ventilation by watching bilateral rise and fall of the chest.				
Retest Approved By: _____			Retest Evaluation: _____		

Evaluator Comments: _____

Candidate Comments: _____

 Evaluator Date

 Candidate Date

 Retest Evaluator Date

 Retest Candidate Date